



Grande Fleur

Grande Fleur cycling a 10-litres glass jacketed reactor

Requirement

This Case Study demonstrates the minimum achievable process temperature and the process temperature control capabilities of the Unistat Grande Fleur when it is connected to an Asahi 10-litre vacuum insulated reactor.

Method

The 10-litres vacuum insulated reactor was connected to Grande Fleur using 1.5-meter M24 metal insulated hoses. The thermofluid used in the system was "DW-Therm". Process control was carried out via a Pt100 sensor located in the process mass. Stirrer speed was set to 100 rpm.

Setup details

Temperature range: -40°C...+200°C 0.6 kW @ +20°C Cooling power:

0.6 kW @ 0°C 0.35 kW @ -20°C

Heating power: 1.5 kW

Hoses: 2*1.5 m metal insulated

HTF: DW-Therm Reactor: Asahi 10-litres glass


jacketed

Reactor content: 7 I DW-Therm Stirrer speed: 100 rpm Control: process Amb. temperature: +25°C

Results

1. Lowest achievable temperature (Tmin):


The graphic shows that a minimum achievable process temperature was -27.8°C.





2. Stability:

The graphic shows the jacket temperature continually adjusting to keep the process temperature constantly at +100°C.



3. Performance:

The table and graphic data show the speed, accuracy and stability of the Grande Fleur as each new set point is reached.

Start T	End T	Approximate time	Av. Ramp Rate	Fastest Ramp Rate
+20°C	-20°C	73 minutes	0.55 K/min	(+10°C to 0°C) 0.7 K/min
-20°C	+100°C	67 minutes	1.8 K/min	(+30°C to +60°C) 2.3 K/min
+100°C	+20°C	71 minutes	1.1 K/min	(+60°C to +30°C) 1.1 K/min

