



Unistat® 925w

Temperature control of the 100l glass jacketed reactor from Chemglass

Requirement

This case study demonstrates the ability of the Unistat 925w to control a 100l glass jacketed reactor from the company Chemglass.

Method

The Unistat and reactor were connected using metal hoses M30. The reactor was filled with M80.100/250.06. Stirrer speed was set to 100 rpm.

Setup details


Temperature range: -90 ... +250°C
 Cooling power: 16,0 kW @ 0°C
 16,0 kW @ -20°C
 15,0 kW @ -40°C
 Heating power: 12,0 kW
 Hoses: 1 x M30 x metal Insulated
 HTF: M90.055/170.02
 Reactor: 100l glass jacketed
 Reactor content: 70l M80.100/250.06
 Reactor stirrer speed: 100 rpm
 Control: Process
 Amb. temperature: +23°C

Results:

1. Lowest achievable temperature in the reaction mass:

The graphic below shows that the minimum achievable temperature of the process mass was -84.5°C. To bring the process mass from -84.5°C back to +20°C took 76 minutes at a rate of 1.4 K/min.

Start T	End T	Time taken	Av. Ramp Rate
+20°C	-84.5°C	228 minutes	0.5 K/Min
-84.5°C	+20°C	76 minutes	1.4 K/min




2. Performance:

Temperature Control from +20°C to -50°C and back to +20°C

The graphics below demonstrate the speed, accuracy and stability of control as the Unistat 925w cycles the process mass between +20°C and -50°C.

Start T	End T	Time taken	Av. Ramp Rate
+20°C	-50°C	64 minutes	1.1 K/Min
-50°C	+20°C	51 minutes	1.4 K/min



Temperature Control from +20°C to -80°C

The graphics below demonstrate the speed, accuracy and stability of control as the Unistat 925w cycles the process mass between 20°C and -80°C.

Start T	End T	Time taken	Av. Ramp Rate
+20°C	-80°C	166 minutes	0.5 K/Min

