

# Unistat® 405w

Cooling a DDPS 2-litre jacketed glass reactor to -20 °C.

### Requirement

This case study looks at the performance of a Unistat 405w cooling a 2-litre glass reactor from 20 °C to -20 °C (40 K) under "process" control.

### Method

The Unistat 405w was connected to the reactor using two 1-metre insulated metal hoses. The reactor was filled with 1.5 litre of "M90.055.03", a silicon based HTF.

### Results

The ramp rate over the temperature change is almost linear at an average speed > 1.7 K/min. taking 23 minutes to reach -20 °C.

### Setup details

Unistat® 405w & DDPS reactor

- Temperature range: -45...250 °C
- Cooling power: 1.3 kW @ 250...0 °C  
0.7 kW @ -20 °C
- Heating power: 1.5 kW / 3 kW
- Pump speed: 3300 rpm
- Hoses: 2x1 m; M24x1.5 (#9325)
- HTF: DW-Therm (#6479)
- Reactor: 2-litre jacketed glass reactor
- Reactor contents: 1.5 litre M90.055.03 (#6259)
- Reactor stirrer speed: 115 rpm
- Control: process

