

## Unistat® 705w

**“Stepping” a Radleys 2-litre “Reactor-Ready” from 20 °C to 250 °C**

### Requirement

Unistats operate over a wide temperature range without the necessity to change the HTF and as they are hydraulically sealed there are no vapours, no oils residue, odour or degradation of the HTF due to oxidation at high temperatures.

This case study demonstrates the ability of a Unistat 705 to heat the process through a series of programmed steps from 20 °C to 250 °C and then cool it back to 20 °C.

### Method

The reactor was filled with 1.5 litre of “P20.275.50” as a thermal load, the stirrer set to 250 rpm and the control to “Process”. The results were recorded on the Huber “Spyware” software.

### Results

The repeated and precise control at high temperatures can be seen in the graphic. Micro-electronics, refrigeration and heating working in unison to provide tight control at elevated temperatures is a practical and here, demonstrated feature of the Unistats wide ranges.



### Setup details

Unistat® 705w & Radleys 2-litre “Reactor-Ready” glass reactor

Temperature range: -70...250 °C  
 Cooling power: 0.6 kW @ -20 °C  
 0.6 kW @ -40 °C  
 0.3 kW @ -60 °C

Heating power: 1.5 kW / 3 kW  
 Pump speed: 3500 rpm  
 Hoses: 2x1 m; M24x1.5 (#9325)  
 HTF: SilOil P20.275.50 (#6157)

Reactor: 2-litre vacuum jacketed glass reactor

Reactor contents: 1.5 litre P20.275.50 (#6157)

Reactor stirrer speed: 250 rpm  
 Control: process

