

## Unistat® 510

**Temperature control of the 100 liters Chemglass reactor**

### Requirement


This Case Study demonstrates the process temperature control and the minimum achievable process temperature when Unistat 510 controls the temperature of the reaction mass in a 100 liters Chemglass reactor.

### Method

The Unistat and reactor were connected using two metal hoses M30. The reactor was filled with 80 liters of DW-Therm. "Process" control was carried out via a Pt100 sensor located in the process mass. Stirrer speed was set to 80 rpm.

### Setup details

Temperature range:	-50 ... +250°C
Cooling power:	5,3 kW @ 0°C 2,8 kW @ -20°C 0,9 kW @ -40°C
Heating power:	6,0 kW
Hoses:	2 x M30 metal Insulated
HTF:	DW-Therm
Reactor:	100 litres glass jacketed
Reactor content:	80 litres DW-Therm
Reactor stirrer speed:	80 rpm
Control:	Process
Amb. temperature:	+29°C




## Results

### 1. Lowest achievable temperature (Tmin):

\*\*In this Case Study, the ambient temperature was a very high 29°C.\*\*

The graphic shows that a minimum process temperature of -38.6°C was achieved.





## 2. Performance:

\*\*In this Case Study, the ambient temperature was a very high 29°C.\*\*

The graphic shows the speed, accuracy and stability of the Unistat 510 as each new set-point is reached.

Start T	End T	Approximate time	Av. Ramp Rate	Fastest Ramp Rate
+20°C	+100°C	114 minutes	0.7 K/min	(30°C to 60°C) 1.4 K/Min
+100°C	+20°C	75 minutes	1.1 K/min	(60°C to 30°C) 1.3 K/Min
+20°C	+150°C	200 minutes	0.7 K/min	(30°C to 60°C) 1.2 K/Min
+150°C	+20°C	107 minutes	1.2 K/min	(+60°C to +30°C) 1.3 K/Min



## 3. Stability:

This Case Study was carried out to simulate more realistic conditions with the Unistat 510 and reactor in full sunlight with an ambient temperature of +29°C.

The graphics show the stability of control at 20°C.

