



Unistat 530w

Unistat 530w controls the simulated exothermic reaction in an uninsulated 100l glass reactor from Chemglass

Requirement

This case study demonstrates the ability of the Unistat 530w to control the process temperature of the reaction mass when a simulated exothermic reaction is taking place in an uninsulated 100l glass reactor from Chemglass.

Method

To simulate the reactions, a 1kW immersion heater was placed into the reaction mass (DW-Therm, 15I) and connected to a controller. The agitator speed was set to 100 rpm. At set points of +70°C and -20°C, the heater was turned "On" at values of 1kW, 750w, 500w and 250w.

Once the temperature had stabilised with the additional heat load of the simulated reaction, the heater was turned "Off". Again, once the temperature had stabilised, the heater was turned "On" again at the new wattage level and the procedure repeated. The results were recorded using a USB thumb drive via the Pilot ONE controller and can be seen below.

Setup details

Temperature range: -55...+250°C Heating power: 12 kW

 Hoses:
 2 x M30 Metal Insulated

 HTF:
 M60.115/200.05

 Reactor:
 glass jacketed reactor 100l


 Reactor content:
 75l M40.165/220.10

Control: process
Stirrer speed: 100 rpm
Amb. temperature: +23°C

Results

1. Simulated reaction at -20°C

| kW | kcal/hr | Process rise | Recovery | Process fall | Recovery | Max Delta T |
|------|---------|--------------|----------|--------------|----------|-------------|
| 1.0 | 860 | 4.3°C | 60 min | 7°C | 32 min | 29K |
| 0.75 | 645 | 3°C | 46 min | 3°C | 26 min | 29K |
| 0.5 | 430 | 2.7°C | 31 min | 2°C | 30 min | 26K |
| 0.25 | 215 | 0.5°C | 32 min | 0.7°C | 32 min | 9K |





2. Simulated reaction at $+70^{\circ}$ C

| kW | kcal/hr | Process rise | Recovery | Process fall | Recovery | Max Delta T |
|------|---------|--------------|----------|--------------|----------|-------------|
| 1.0 | 860 | 1.6°C | 26 min | 1.6°C | 27 min | 15.5K |
| 0.75 | 645 | 1.5°C | 27 min | 1.6°C | 25 min | 14.8K |
| 0.5 | 430 | 1.3°C | 27 min | 1.4°C | 28 min | 13.1K |
| 0.25 | 215 | 0.9°C | 25 min | 1°C | 27 min | 8.1K |

