



Setup details

Unistat® 825w & Buchi Glas Uster «miniPilot» 10 reactor

- Temperature range: -85...250 °C
- Cooling power: 2.4 kW @ 0...-40 °C
1.5 kW @ -60 °C
- Heating power: 3.0 kW
- Pump speed: 3500 rpm
- Hoses: 2x1.5 m; M30x1.5 (#6386)
- HTF: DW-Therm (#6479)
- Reactor: 10-litre jacketed glass reactor
- Reactor contents: 7.5 litre M90.055.03 (#6259)
- Reactor stirrer speed: 400 rpm
- Control: process

Unistat® 825w

Heating and cooling a Buchi Glas Uster 10-litre reactor

Requirement

The graphic illustrates the performance temperature of a Unistat 825w heating and cooling a Buchi Glas Uster 10-litre reactor between 20 °C and 60 °C.

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 7.5 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The process temperature heats at a ramp rate of 3.8 K/min. taking 30 minutes to reach the set-point. The cooling ramp rate is at a rate of 5.7 K/min. and takes 27 minutes to cool through 40 K (60 °C to 20 °C).

