CC-106A

Clear View Circulator Bath consisting of transparent baths made of polycarbonate with temperature range up to max. 100°C. Powerful pressure and suction pump made of high-resistant plastic. With adjustable overtemperature protection according to DIN 12876.

Pilot ONF:

The new Pilot ONE controller with pioneering technology and advanced control functions brings numerous advantages to routine work. The extensive features list includes a brilliant 5,7" TFT touchscreen display, USB and network connections, an integrated technical glossary and language support in 13 languages (EN, DE, FR, IT, ES, RU, CN, PT, JP, CZ, PL, KO, TR). The Pilot ONE has a convenient navigation system with easily remembered icons and menu categories which are colour sorted to make routine work simpler. Thanks to a favourites menu and One-Click operator guidance all important information is always just a few keystrokes away. Software wizards also help you to set up, ensuring correct settings. The USB port allows connection of the system to a PC or notebook. Together with the Spy software, requirements such as remote control or data transmission are easily achieved in a cost-effective manner. Network integration is easy with the internet port.

The range of functions can be expanded very easily via E-grade at any time by entering a unit specific upgrade code:

E-grade "Exclusive": TAC (True Adaptive Control) - self optimising internal and cascade control, selectable temperature control mode (Internal/Process), programmer with 3 programs (max. 15 steps), ramp function (linear), 5 point calibration, scalable graphic display, favourites menu, display resolution 0,01 K.

E-grade "Professional": Programmer with 10 programs (max. 100 steps), ramp function for temperature gradients (linear and non-linear), 2nd set point, user menus (Administrator level), calendar start.

4-year warranty - registration required.

Technical data according to DIN 12876

Operating temperature range with water cooling with refrigerator Temperature stability at 70°C temperature set point / display Absolute accuracy Internal temperature sensor Sensor external connection

Safety classification Heating power at 240V Heating power at 230V Heating power at 220V Heating power at 208V Heating power at 200V Pressure pump

Interface digital

max. delivery max. delivery pressure Suction pump max. delivery (suction)

max. delivery pressure (suction) Pump connenction (optional)

Bath volume min. filling capacity

Width bath opening WxD/ bath depth Overall dimensions WxDxH **

Net weight

Power supply requirement

max. current min. Fuse max. Fuse Pressure equipment category Degree of Protection min. ambient temperature max. ambient temperature

Peter Huber Kältemaschinenbau SF

25...100 °C 20...100 °C

15...100 °C 0.02 K

5,7" colour Touchscreen setup for calibration

Pt100 Pt100

Ethernet, USB (Host u.

Device), RS232

III / FL 2.1 kW 2 kW 1,8 kW 1,6 kW 1,5 kW

27 I/min 0.7 bar yes 22 I/min 0.4 bar M16x1 male

4.4 [251

130x110/ 150 mm 147x307x330 mm

5 kg

200-240V 1~/2~ 50/60Hz

10A 16A Art. 4.3 PED

IP20 5°C

40 °C

Werner-von-Siemens-Str 1

Tel 0781/9603-0 Fax 0781/57211 www huber-online com

Order-No.: 2049.0001.01

10 A

D-77656 Offenburg

from Serial-No.: 395403 1.0/20

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions. Illustrations can deviate from the original.

Accessories and periphery: mini-USB cable #54949*, Bath cover #37533, Bath Bridges #19592, Cooling coil #30554, pump adaptor #19606, Note: When using Huber pump adapter: Polyglycol is not permissible to be used as a heat transfer fluid,polycarbonate test tube racks Typ A-F, nozzle #33288, DS level regulator #9580

* standard equipment

Output data valid for: Room temperature 20°C

In accordance with EN60034-1 the following voltage and frequency tolerances are valid:

Voltage + / - 10%, as long as the frequency tolerance does not run in the opposite direction.

Example: -10% voltage and +3% frequency -> not allowed!

-10% voltage and -3% frequency -> allowed.

Information to Electromagnetic compatibility:

Classification (disturbance) to EN55011: Class A, Group 1

Standard delivery conditions - Power cable configuration:

- 1. Single / two-phase devices (100V to 240V) --> with power cable and country-specific plug (please specify when ordering)
- 2. Three-phase devices with current consumption less than 63A --> with cable, without plug
- 3. Three-phase devices with current consumption greater than 63A --> without cable, without plug
- ** Please respect space requirements. See operating conditions at www.huber-online.com

Peter Huber Kältemaschinenbau SE Werner-von-Siemens-Str. 1 D-77656 Offenburg Tel 0781/9603-0 Fax 0781/57211 www.huber-online.com