huber

CC-308B

Heating Circulator Bath with housing, bath and all moistened parts are made of stainless steel. With cooling coil for water-cooling (3/8"), pressure- and suction pump made of high-resistant plastic. With adjustable overtemperature protection according to DIN 12876.

Pilot ONF:

The new Pilot ONE controller with pioneering technology and advanced control functions brings numerous advantages to routine work. The extensive features list includes a brilliant 5,7" TFT touchscreen display, USB and network connections, an integrated technical glossary and language support in 13 languages (EN, DE, FR, IT, ES, RU, CN, PT, JP, CZ, PL, KO, TR). The Pilot ONE has a convenient navigation system with easily remembered icons and menu categories which are colour sorted to make routine work simpler. Thanks to a favourites menu and One-Click operator guidance all important information is always just a few keystrokes away. Software wizards also help you to set up, ensuring correct settings. The USB port allows connection of the system to a PC or notebook. Together with the Spy software, requirements such as remote control or data transmission are easily achieved in a cost-effective manner. Network integration is easy with the internet port.

The range of functions can be expanded very easily via E-grade at any time by entering a unit specific upgrade code:

E-grade "Exclusive": TAC (True Adaptive Control) - self optimising internal and cascade control, selectable temperature control mode (Internal/Process), programmer with 3 programs (max. 15 steps), ramp function (linear), 5 point calibration, scalable graphic display, favourites menu, display resolution 0,01 K.

E-grade "Professional": Programmer with 10 programs (max. 100 steps), ramp function for temperature gradients (linear and non-linear), 2nd set point, user menus (Administrator level), calendar start.

4-year warranty - registration required.

Technical data according to DIN 12876

Operating temperature range with water cooling

with refrigerator

Temperature stability at 70°C temperature set point / display

Absolute accuracy

Internal temperature sensor

Sensor external connection

Interface digital

Safety classification Heating power at 115V Heating power at 100V

Pressure pump max. delivery

Delivery pressure (head)

Suction pump

max. delivery (suction)

Delivery suction pressure (head)

Pump connection

max. permissible kin. viscosity

min. filling capacity Bath volume

Bath capacity with displacement rack

Width bath opening WxD

Bath depth

Height of bath opening
Overall dimensions WxDxH **

Net weight

Power supply requirement

max. current

Pressure equipment category

Degree of Protection min. ambient temperature max. ambient temperature 28...300 °C

20...300 °C -20...300 °C

0.02 K

5,7" colour Touchscreen

setup for calibration

Pt100 Pt100

Ethernet, USB (Host u.

Device), RS232

III / FL 1,5 kW 1 kW yes 25 l/min 0,7 bar

yes 18,5 l/min 0,4 bar M16x1 male

50 mm²/s 5,8 I 7,6 I 5.2 I

130x110 mm 155 mm 194 mm

242x404x392 mm

18 kg

40 °C

100-115V 1~ 50/60Hz

14 A 15 A Art. 4.3 PED IP20 5 °C



Order-No.: 2006.0003.01

Fax 0781/57211

www huber-online com

from Serial-No.: 1.2/20

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions. Illustrations can deviate from the original.

Included Accessories:

mini-USB cable #54949, hose connector NW12, blank plug, sleeve nuts thread,

Optional accessories:

connecting hoses, drain valve , displacement insert, temperature control / - connection hoses, thermofluids, further accessories, etc.: see catalog.

Output data valid for: Room temperature 20°C

In accordance with EN60034-1 the following voltage and frequency tolerances are valid:

Voltage + / - 10%, as long as the frequency tolerance does not run in the opposite direction.

Example: -10% voltage and +3% frequency -> not allowed!

-10% voltage and -3% frequency -> allowed.

Information to Electromagnetic compatibility:

Classification (disturbance) to EN55011: Class A, Group 1

Standard delivery conditions - Power cable configuration:

- 1. Single / two-phase devices (100V to 240V) --> with power cable and country-specific plug (please specify when ordering)
- 2. Three-phase devices with current consumption less than 63A --> with cable, without plug
- 3. Three-phase devices with current consumption greater than 63A --> without cable, without plug
- ** Please respect space requirements. See operating conditions at www.huber-online.com

Peter Huber Kältemaschinenbau SE Werner-von-Siemens-Str. 1 D-77656 Offenburg Tel 0781/9603-0 Fax 0781/57211 www.huber-online.com