CC-902





Refrigerated Heating Circulator Bath with air-cooled cooling machine. Powerful, variable speed, pressure and suction pump, evaporator (cooler) and housing of stainless steel, CFC and H-CFC free. With adjustable overtemperature protection according to DIN 12876.

Pilot ONF:

The new Pilot ONE controller with pioneering technology and advanced control functions brings numerous advantages to routine work. The extensive features list includes a brilliant 5,7" TFT touchscreen display, USB and network connections, an integrated technical glossary and language support in 13 languages (EN, DE, FR, IT, ES, RU, ZH, PT, JA, CS, PL, KO, TR). The Pilot ONE has a convenient navigation system with easily remembered icons and menu categories which are colour sorted to make routine work simpler. Thanks to a favourites menu and One-Click operator guidance all important information is always just a few keystrokes away. Software wizards also help you to set up, ensuring correct settings. The USB port allows connection of the system to a PC or notebook. Together with the Spy software, requirements such as remote control or data transmission are easily achieved in a cost-effective manner. Network integration is easy with the internet port.

The range of functions can be expanded very easily via E-grade at any time by entering a unit specific upgrade code:

E-grade "Exclusive": TAC (True Adaptive Control) - self optimising internal and cascade control, selectable temperature control mode (Internal/Process), programmer with 3 programs (max. 15 steps), ramp function (linear), 5 point calibration, scalable graphic display, favourites menu, display resolution 0,01 K.

E-grade "Professional": Programmer with 10 programs (max. 100 steps), ramp function for temperature gradients (linear and non-linear), 2nd set point, user menus (Administrator level), calendar start.

4-year warranty - registration required.

The dimensions specified below are approximate and may be subject to change.

Technical data according to DIN 12876

-90200 °C	
0,02 K	
5,7" colour Touchscreen	Order-No.: 2026.0015.01
Pt100	
Pt100	
Ethernet, USB (Host u. Device), RS232	
III / FL	
1,2 kW	
1,2 kW	
1,2 kW	
1,2 kW	
1,1 kW	
0,9 kW	
0,6 kW	
0,2 kW	
air-cooled, natural refrigerant	
R-1270 (A3, H220)	
0	
R-170 (A3, H280)	
0,437	
UN 3358	
yes	
25 l/min	
0,7 bar	
yes	
18,5 l/min	
0,4 bar	
M16x1 male	
50 mm²/s	
5 I	
41	
	0,02 K 5,7" colour Touchscreen Pt100 Pt100 Ethernet, USB (Host u. Device), RS232 III / FL 1,2 kW 1,2 kW 1,2 kW 1,2 kW 1,1 kW 0,9 kW 0,6 kW 0,2 kW air-cooled, natural refrigerant R-1270 (A3, H220) 0 R-170 (A3, H280) 0,437 UN 3358 yes 25 I/min 0,7 bar yes 18,5 I/min 0,4 bar M16x1 male 50 mm²/s 5 I

Technical data according to DIN 12876

Width bath opening WxD / bath depth 120 x 110 / 150 mm

Height of bath opening 712 mm

Overall dimensions WxDxH **

Power supply requirement

Degree of Protection

max. ambient temperature

min. ambient temperature

550x600x911 mm

208V 2~ 60Hz

IP20

40 °C

5 °C

from Serial-No.: 1.4/21

Technical details and dimensions are subject to change. No liability is accepted for errors or omissions. Illustrations can deviate from the original.

Output data valid for: Room temperature 20°C. If the ambient temperature rises, the cooling capacity may drop.

in accordance with EN60034-1 the following voltage and frequency tolerances are valid:

Voltage + / - 5% with a simultaneous frequency tolerance of + / - 2%

Example -5% voltage and + 2% frequency -> not allowed!

-5% voltage and - 2% frequency -> allowed

Information to Electromagnetic compatibility:

Classification (disturbance) to EN55011: Class A, Group 1

Standard delivery conditions - Power cable configuration:

- 1. Single / two-phase devices (100V to 240V) --> with power cable and country-specific plug (please specify when ordering)
- 2. Three-phase devices with current consumption less than 63A --> with cable, without plug
- 3. Three-phase devices with current consumption greater than 63A --> without cable, without plug

This equipment is compliant to US-SNAP and all applicable EU laws. The US-SNAP end-use for this equipment is the industrial process refrigeration. Certification by a Notified Body upon request.

** Please respect space requirements. See operating conditions at www.huber-online.com

Peter Huber Kältemaschinenbau SE Werner-von-Siemens-Str. 1 D-77656 Offenburg Tel 0781/9603-0 Fax 0781/57211 www.huber-online.com